

PECTUS BAR SERIESE


PECTUS BAR FOR EXCAVATUM/
CARINATUM DEFORMITY



PECTUS BAR SERIES

Product	Part No.	Purpose
Pectus bar	TIS1031	Funnel chest, Pectus excavatum, Costi-cartilage correction
Irrotatable bar	TIS1031 N TIS1031 P	
Pigeon bar	TIE1131	Pectus carinatum


TIS 1031 PECTUS BAR – SINGLE BAR



TIS 1031 N(ABOVE) / TIS 1031 P(BELOW) – DOUBLE BAR





TIE 1131 PIGEON BAR




PECTUS BAR SERIES



TIS 1031 PECTUS BAR (SINGLE BAR)



TIS 1031 N(ABOVE) / TIS 1031 P(BELOW)



TIE 1131 PIGEON BAR



PECTUS BAR SERIES

PECTUS BAR FOR EXCAVATUM/CARINATUM DEFORMITY



TIS 1031 SINGLE BAR



TIS 1031 DOUBLE BAR



TIE 1131 PIGEON BAR

인증된
티타늄소재
:Ti6Al4V ELI

DESIGNED
INSTRUMENTS
FOR
PRECISENESS
AND EFFICIENCY
OF SURGERY

INNOVATIVE
STABILIZER

: EASY INSERTION
BAR INTO ANGLED
POSITION

VARIOUS
LENGTHS

: TIE1131 (12 TYPE)
TIS1031 (12 TYPE)
TID10AI (ABOVE: 23,
BOTTOM:12)

DOWNSIZED
STABILIZER

: REDUCE INSERTION
TIME AND SOFT
TISSUE IRRITATION

ROUNDED
AND BLUNT
EDGES OF BAR

: MINIMIZE SOFT
TISSUE DAMAGE



Sum up for PECTUS BAR character

- ① Titanium bar – light, high tension, safe(no allergy)
- ② round edge
- ③ TID - No need stabilizer, better cosmetic effect

PECTUS BAR FOR EXCAVATUM/CARINATUM DEFORMITY



TIS 1031 PECTUS BAR



TIS 1031 N(ABOVE)
TIS 1031 P(BELOW)



TIE 1131 PIGEON BAR

